

Varilux® Elipse™

80/85mm
75/80mm
70/75mm





Fig. A. Fitting heights positioned too low can result in loss of useable reading area - one of the most common causes of headaches. To assure proper postoperative fitting, cross marks over the horizontal. Check to be sure that the vertical fitting height is a minimum of 1.4 mm.

© 2004 Essilor of America, Inc.
Varilux is a registered trademark of Essilor International, S.A.
Elipse is a trademark of Essilor of America, Inc.



Minimum suggested fitting height for Elipse is 1.4 mm.

LVAR 30018
0804CP