


# HOYALUX SUMMIT ECP LAYOUT GUIDE



The addition is engraved on the temporal side and the design/material code is engraved on nasal side.

Lens Design/Material Code

SP 1	=	1.50 CR-39
SPX	=	1.53 Phoenix
SP 3	=	1.60 EYAS
SP 4	=	1.70 EYRY


HOYALUX  
**summit** ecp

The EYRY 1.7 Index Summit ECP is designed with 2.5mm decentration


**HOYA**

# Datum Line Chart

L



R



# HOYA

MILLET 03/06 377142